

Automated External Defibrillator Inspection Log



Company	
Date	
Make & Model	
Serial Number	
Inspector	
Location	

Criteria	Yes	No	Corrective Action / Comments
Is the status indicator light green or flashing green?			
Is the AED "chirping or beeping"?			
Are the pads within the expiration date?			
Is the battery within the expiration date?			
Is the AED visibly damaged?			
Is the AED clear of debris and moisture?			
Is a first aid kit available with a razor, gloves, CPR mask, antiseptic towelettes, & scissors?			
Are emergency postings present and visible?			
If there is an AED wall cabinet with an alarm, is the door alarm working? (optional)			

*AED Log makes it easy to track AED inspections, expiration dates, and program compliance online.
Sign up at aedlog.com to manage 1 AED for free, or upgrade for unlimited AED management.*